- To register program **DSA Monitor**:
- 1. Start the application.
- 2.Use menu command: Help \rightarrow Registration:



to open registration window:

🜌 DSA Monitor Registration		×
Program was registered by:		
SONOPAN Sp. z o.o. 15-950 Białystok ul. Ciołkowskiego 2/2		
Features:	VEC	
Data visualization in tables: Data visualization on charts:	YES YES	
Statistical calculations:	YES	
 Data creation by user: Import data from other SONOPAN files: 	YES YES	
Data import from any compatible file:	YES	
Measured data edition: Full data access:	YES YES	
🔗 🖉 Register 🕐 About reg	jistration	X Close

3. Button Register: Pegister will open dialog, in which registration key file DSA Monitor. RegKey have to be selected.